

ULTRA LIGHT OPTICS™

The **SMALLEST** and **LIGHTEST** LED Headlight

"Experience superior comfort & illumination"



Battery Pack Specifications

Dimensions: 3.0" x 1.0" x 2.5"

Weight: 5oz

Battery: High Capacity Lithium Ion

Charge time: approx 3 hours

Battery Life: approx 8 hours
est 400-700 charges

Type: Regulated, steady power

1 year limited warranty on the battery



SPECIAL OFFER
RM 3,500

Light Specifications

Diameter: 0.65" (smaller than a dime)

Length of Module: 0.76"

Weight: approx 1.5g, or 0.05oz (lighter than a nickel)

Intensity: approx 3,600 footcandles (39,000 Lux)

Light Color: 6,000k White Light

Continuous run time: approx 8 hours
consistent brightness

Spot Size: 2.5" @ 12"

Longevity: up to 50,000 hours

*All prices are subject to 6% GST

*Terms & Conditions Apply



Customized mounting brackets

Orasoptic™

(Flip-Up) with Curvy Bar

(Flip-Up) with Flat Bar

Any (TTL)

Surgitel™

Any (Flip-Up) Surgitel (TTL)

Flak Jacket (TTL) Radar (TTL)

Half Jacket (TTL) A-Wire (TTL)

Aero (TTL) M-Frame (TTL)

Design for Vision™

Nike Curfew (TTL)

Small Sport Frame (TTL)

Nike Pounce/Revolve (TTL)

Black/ Buddy Holly/Yeoman (TTL)

Titanium/Chocolate (TTL)

Nose Bridge Look: = =

Nike Skylon Ace (TTL)

Safety (TTL)

Sky View

Heine™

Brasseler™

Sheer Vision™

Carl Zeiss™

Perioptix™

Q-Optics™

Synca™

Shield™

(Flip-Up), (TTL)

Any Flip-Up, Under Armour™ (TTL)

(Flip-Up Prismatic) (Flip-Up Galilean)

(Flip-Up) with Curvy Bar

(Flip-Up) with Flat Bar

Adidas Adivista (TTL)

(Flip-Up), (TTL), Uvex (TTL)

Exam Vision

Clip-On | Drill-In

OVI (Flip-Up)

Medlite (Flip-Up)

Rose (Flip-Up)

Keeler (Flip-Up)

Eagle (Flip-Up)

Sandy Grendel

Reading Glasses

Safety Goggles

Headband

All trademarks are the properties of their respective owners



TITAN DENTAL SOLUTIONS SDN BHD (1211585-W)

Tel: 03-6206 1175 Sales: 012 731 1529

No. 20-2, Jalan 27/70A, Desa Sri Hartamas 50480 Kuala Lumpur, Malaysia.

TitanDentalSolutions

Please call to order your replacement parts or accessories

 TitanDentalSolutions



Standard Battery



Advanced Sensor Battery



Advanced Touch Medical Battery



Cordless Hybrid



Micro USB Charger



Universal Charger



Wire Extension



ELITE Wire Replacement



Composite Filter








)



Miscellaneous Parts

Customized Adaptors

Available Colors

-  Silver
-  Black
-  Pink
-  Red
-  Royal Blue



Sport Frame Goggles



Premier Case



Standard Case



Deluxe Case